



High Containment Isolators

Bringing an advanced level of containment to your work environment



Howorth UltraSAFE High Isolators are designed utilising our specialist experience including stringent safety measures to ensure protection of personnel and the environment.

- Designed to safely handle potent compounds from OEB3 to OEB6
- Negative pressure cascade
- Single chamber design
- Material Pass Box, RTPs and Bag Out Systems
- Full ergonomic assessment
- Operational flexibility
- Designed for safe or hazardous area classification
- Safety features designed into the system including interlocks, purge systems and safe change systems
- CIP/WIP Capabilities
- PLC control system including interface with site networks.
- Compliance with all appropriate design standards and guidelines



Operational flexibility

Work area for installation of various laboratory or weighing equipment
Integrated Power Sockets and Service outlets
Off floor storage



Interlocked Pass Box

Interlocked pass box to allow safe transfer of materials and tools during process operations



Wash in Place (WIP)

Hand held wash hose providing manual cleaning of the chamber internals
Spray balls are also available if required.

Head Office

Howorth Air Technology Limited
F2C Multiply, Lomax Way,
Logistics North, Bolton,
BL5 1FQ,

Contact Us

info@howorthgroup.com
+44 (0) 1204 571 131

USA Office

Howorth Air Technology Inc.
4601 Corporate Drive NW, #140
Concord NC 28027
usa@howorthgroup.com